
LabActor

Main control panel



[abc] Monitor

List of active measurement modules

III Start

Start measurement

TFI Stop

Stop measurement

TFI View

View Monitor data file

TFI Read cfg

Read application configuration from file

TFI Write cfg

Write application configuration to file

III Monitor configuration

Open Monitor configuration

TFI Exit

Exit application

TFI Clear log

Clear error log

Monitor Running

Measurement state


Monitor File

Latest Monitor data file

Error log

Monitor

Monitor configuration



Configuration

Monitor configuration

DBLI Update period

Request data collection this often

1321 Repeat times

Complete this many measurements then stop (or keep going if <0)

TFI Open plot


Open new Monitor data plot on start

IDBL! File defragment period

Defragment Monitor data file this often when running

Spectrum Recorder

Spectrum Recorder module sweeps frequency and records complex response using Red Pitaya



Essile Configuration

Spectrum recorder configuration

DBL Start frequency

Spectrum starting frequency

BL End frequency

Spectrum ending frequency

DBLI Frequency step

Minimal frequency increment

DBL Point interval

Interval between frequency increments

Address

Network address of Red Pitaya

TFI Open plot

Open new spectrum plot window at the start of measurements

[DBL] Focus region centers

If not empty, increase frequency steps outside listed regions

DBL Regions width

Width of focus regions

Out-of-focus step acceleration

If > 1, increase frequency steps outside listed regions by this factor

TFI View XY Plot

View latest spectrum file

Spectrum plot

Latest spectrum

Spectrum File


Latest spectrum file

Sweep duration

Estimated duration of single spectrum recording, sec

Line Detector

Line Detector module finds amplitude peaks in spectrum file and fits complex Lorentzian around peaks



Input file

Spectrum file to process (auto-updates when monitor is running)

ESSIP Configuration

Line Detector configuration

Initial fit guess

Run peak finder to guess fit parameter or use Best fit results values

IDBL Fit window

Width of complex Lorentzian fit window around the peak

BLI Re-sample

If > 0, re-sample signal to this frequency step

DBL Peak width guess

Guess width for peak finder

Finder algorithm

Peak finder type

Background removal

Background removal type for peak finder

Smoothing

Smoothing type for peak finder

DBLI Smooth window

Smoothing amount for peak finder

DBL Relative threshold

Smallest peak magnitude relative to spectrum maximum

Belling Bgd window

Background removal width for peak finder

Spectrum recorder name

Name of the spectrum recorder module to watch

Follow peaks

Auto-update Spectrum recorder focus regions with best fit results

TEI Create rusin.dat

Auto-create rusin.dat

шы ...with header

Prepend rusin.dat with this header

Re-process spectrum selected in the Input file window

Spectrum plot

Spectrum data with complex Lorentzian fit

Peak Plot


Filtered and smoothed magnitude with peaks found

[DBL] Best fit results

Latest complex Lorentzian fit results. Manual edits allowed when used with "Last fit" guess mode

Data Plot Module

Plots listed TDMS data files



- File(s)
 - List of files to plot. Click browse buttons on the right to add/replace
- X column

Selects X-axis channel name. Avaliable names determined by top file on the lsit.

Y column(s)

Selects Y-axis channel name(s). Hold Ctrl or Shift for multi-channel selection. Avaliable names determined by top file on the lsit.

Plot type

Switch between X-Y and Waveform style plots

File history depth

Keep up to this many previous files in the list when new measurement file arrives